


FOR NON-SURGICAL
BODY SHAPING
and SKIN PERFECTION



▶ **ENDOSPHERES THERAPY®** is a brand-new treatment harnessing the power of **COMPRESSIVE MICROVIBRATION®**, patented system engineered for non-invasive skin rejuvenation, body shaping and physical therapy. The method has been studied by **influential medical research centers**, that have certified the efficacy of the treatment.

▶ The treatment is delivered using a cylinder made up of special **55 silicone spheres** positioned in a honeycomb-like pattern, that have a specific density and diameter.




▶ **"HONEYCOMB"** arrangement

5 ACTIONS


The device uses **COMPRESSIVE MICROVIBRATION** to transmit low frequency vibrations creating a pulsed rhythmic action on the tissue. The wave resulting from the vibrations **has 5 therapeutic benefits on the tissue**. It provides a well-rounded treatment. Initially it restores the normal state of tissues and microcirculation and in later stages works on reshaping and contouring the body.

pain
RELIEF




Compressive Microvibration acts on mechanoreceptors to reduce pain and tissue inflammation. It improves tissue oxygenation, reducing painful swelling.

VASCULAR
action




The honeycomb arrangement of the spheres stimulates the vascular system, affecting the blood flow in the tissues. The **"vascular workout"** of the tissue improves microcirculation.

LYMPHATIC
drainage




Compressive Microvibration uses the rotational movement of the spheres and the direction of the cylinder on the tissue to activate the lymphatic system to work correctly. This helps to **eliminate excess fluids that have accumulated as a result of poor circulation**.

body
CONTOURING




The oscillations of the spheres create a **push pull action** which remodels the tissue. The shape and contour of the body is improved after only a few sessions, as well as reducing unsightly orange peel effect visible on skin.

MUSCLE
toning



Compressive Microvibration technology tones the muscles as they actively resist the compression. The vibration stimulates the muscles and can improve tone in the treated areas.




excellent

▶ **VALUE FOR MONEY**
for a total body device

▶ **VERSATILE TECHNOLOGY**
with safe and reliable results

▶ **NEW GENERATION**
TECHNOLOGY



Low frequencies mechanical oscillations
work in depth on the main causes
of cellulite.

▶ **ENSURES**
customer retention



FACIAL TREATMENT *Areas of application*

Mouth • Eyes • Forehead • Cheeks
Cheekbone • Glabella
Neck and lower neck area • Breast
Pre and post Injectable treatments

- ▶ SMOOTHES WRINKLES
- ▶ IMPROVES VASCULARIZATION AND OXIGENATION
- ▶ LIFTS AND FIRMS THE SKIN
- ▶ RELAXES CONTRACTED MUSCLES (EXPRESSION LINES)
- ▶ STIMULATES COLLAGEN PRODUCTION



▶ **COMFORTABLE**
for patients

▶ **100%**
non-invasive

▶ **VISIBLE RESULTS**
noticed after only
a few sessions



BODY TREATMENT *Areas of application*

Cellulite • Orange Peel skin • Contracted Muscles • Oedema
Swollen legs • Arms • Buttocks • Ankles • Abdominal fat
Breast tightening • Post-pregnancy treatment

- ▶ REDUCES LYMPHATIC FLUID RETENTION
- ▶ REDUCES ORANGE PEEL SKIN
- ▶ IMPROVES VENOUS CIRCULATION
- ▶ IMPROVES BODY SHAPE
- ▶ TONES MUSCLES



▶ **INCREASE OF**
salon/clinic revenue

▶ **EXCELLENT**
return on investment

Clinical observations have confirmed **Endosphères Therapy** as a **valid and trustworthy treatment** providing a solution for cellulite, poor venous circulation, fluid retention and muscle tension.

In the studies, in order to obtain the correct diagnosis, patients were subjected to a clinical evaluation. This included the inspection and palpation of the tissue. Echo colour Doppler, echography, thermography and biometrical impedance analysis were carried out. Analysis was carried out after every six sessions to give an insight into the tissue changes. The results show the efficacy of **Compressive Microvibration**.

Courtesy
"G. d'Annunzio"
UNIVERSITY
CHIETI/PESCARA - ITALY
Dept. Medicine Physical
Therapy and Rehabilitation

CLINICALLY PROVEN RESULTS

Results after **8** sessions



Results after **15** sessions



Results after **18** sessions



Results after **3** sessions




Results after **3** sessions



Results after **3** sessions



Immagini: Archivio XI Congresso Internazionale SIES Medicina e Chirurgia Estetica



EOS

TECHNICAL DATA

Diameter	315mm(up)/364mm(down)
Height	1.150 mm
Weight	16 kg
Vibration	29 - 310 Hz
Cooling	ED/AIRCOOLING
Touch Screen	Capacitivo 7 Colour
Sound Level	56dBA
Electrical Features	da 100 a 240 V./ 50-60Hz

www.endospheres.com

The Endosphères Therapy EOS patent
belongs to the Fenix Group and is manufactured by the Fenix Group



The Endosphères Therapy brand belongs to the Fenix Group, a 100% Italian company. The group has been at the forefront in the field of aesthetics and medical-aesthetic equipment, always basing its strength on Research and development.

The Fenix Group has developed Endosphères Therapy, a new technology, using a patented Compressive Microvibration method. In recent years this has proved itself as the treatment of excellence to combat edematous fibrosclerotic panniculopathy (cellulite), adipose accumulations and face rejuvenation.

The Fenix Group is an innovative company investing in

marketing and PR campaigns both in trade press and consumer press. It is active across all social media platforms. Customers are assisted with in salon/clinic marketing to promote the treatments locally in their area. Service and technical support is available throughout the world.

The Fenix Group's mission is to design high quality electro-medical devices and manufacture the most innovative, effective and science-based technologies. These essential features ensure safety and efficacy making the Fenix Group a successful and customer-driven company.



Headquarters ITALY

Fenix Group s.r.l. - Strada Lungofino, 187 - c/o C.C. Ibisco B9/10 Città Sant'Angelo (PE)
Tel. **+39 085 4483235** - info@fenixgroup.it



Branch Office U.K.

Fenix Group LTD - 51 Holland Street, London W8 7JB
Tel: **+44(0)2070960040** - ukmanagement@fenixgroup.it



Branch Office SPAIN

Fenix Endospheres Therapy SL - C/ Cirilo Amoros, 6 Entresuelo 46004 Valencia
Tel. **+34 627 508380** - espmanagement@fenixgroup.it



Branch Office FRANCE

Fenix France S.A.S. - Startway 43-47 Avenue de la Grande Armee - 75016 PARIS
Tel: **+33(0)170395146** - frmanagement@fenixgroup.it

For other countries look for distributors on www.endospheres.com